

*We've raised our standard  
with no additional cost to you!*

## **Snap-Lock Expander**



**The leader in quality  
and value.**

 **NorthStar**  
ORTHODONTICS, INC

# NorthStar

ORTHODONTICS, INC

218 Industrial Park Road West  
P.O. Box 146  
Park Rapids, MN 56470

[www.northstardental.com](http://www.northstardental.com)

1-800-346-0011

1. The Snap-Lock Expander does not turn back unless you or your patient turns it back. After the activation of the spindle, the flat spring snaps in.
2. The hole in the spindle is always in the right place for reactivation.
3. When the Snap-Lock Expander is properly activated by 1/4 turn you will be able to feel the screw lock.
4. If overactivated, the Snap-Lock Expander can be turned back.
5. You can control any situation without concerns about the appliance not doing what it is supposed to during treatment.
6. The Snap-Lock Expander will not fall apart when fully expanded and the spring prevents the key from impinging on the palatal tissue.
7. The spindle is oriented symmetrically within the housing and the expansion key hole preventing unwanted poking in the palate.

**\$10 OFF**

one

**Snap-Lock  
Expander**

*with return of this card*



Recycled Paper